

CTAB DNA Extraction Protocol

(3/29/00)

1. Sample leaves and two grinding beads in 2 ml tube.
2. Add 600 μ L CTAB and 12 μ L β -mercaptoethanol in 2mL tube in hood.
(If you have CTAB mixed with β -mercaptoethanol, add only 600 μ L)
3. Grind sample in the grinding machine (3-5 min; 13.3K rpm).
4. Place in 65°C for 30-40 minutes.
5. Add equal volume (600 μ L) 24:1 chloroform / isoami mix to each tube and mix. Solution should be cloudy. Shake the tubes well for mixing.
6. Centrifuge for 10 minutes @13K.
7. Transfer upper phase (400ul) into new 1.5mL tube containing 280 μ L cold isopropanol (found in freezer). Do not transfer any chloroform into new tube. Shake the tubes well for mixing.
8. Place tubes on ice for 5 minutes.
9. Centrifuge for 15 minutes @13K.
10. Dump supernatant.
11. Add 400 μ L cold 70% ethanol (in freezer) to wash DNA pellet.
12. Centrifuge for 5 minutes @ 13K.
13. Decant supernatant and dry the pellet in vacuum dryer 20-30 minutes, until dry.
14. Resuspend DNA pellet in 200 μ L ddH₂O. Quantify with nanodrop and dilute to desire concentration.
 - >If doing PCR same day, leave tubes at room temperature for 1 hour.
 - >Otherwise, store DNA solutions at -20°C /4°C.